# FAQ - Basic Research

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### Microdialysis Technique

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<b>Question:</b>	What is recovery? 什么是回收率？
------------------	------------------------------

**Answer:** The recovery of a particular substance is defined as the concentration in the dialysate expressed as percent of the concentration in the interstitial fluid.

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<b>Question:</b>	How does the length of the membrane affect the recovery? 膜长度对回收率的影响？
------------------	---

**Answer:** A longer membrane and a lower flow rate will give a higher recovery.

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<b>Question:</b>	How shall the perfusion fluid be composed? 灌流液应如何配制？
------------------	---

**Answer:** Ideally it should be as close as possible to the composition of the extracellular fluid. However, you may want to change the concentration of sodium, potassium or calcium in order to influence the membrane function in the region you are studying. CMA can offer Perfusion Fluid T1 for peripheral tissue and Perfusion Fluid CNS for brain tissue.

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**Answer:** The CNS perfusion fluid is intentionally not buffered in order to allow it to take on the same pH as the brain's interstitial fluid. The pH of a non-buffered solution varies between 5-8.

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**Answer:** Since our perfusion fluid is unbuffered it can hardly itself impact the pH in the tissue. The perfusate will take on the same pH as the surrounding tissue. But the buffered substances in the tissue can have different recovery over the membrane which indirectly can have a small effect on the pH of the tissue.

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The ability of molecules to pass the membrane decreases logarithmically with the increase in molecular weight. By experience we know that most substances with a molecular weight up to 5,000 Da can be dialyzed when using a 20,000 Da membrane. This is of course very dependent on the substance and the sensitivity of the analytical method.

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## Probes

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**Question:** Which length of the microdialysis probe membrane shall I use?  
微透析探针的膜的长度有哪些？

**Answer:** A longer membrane gives a better recovery of the substances you are interested in but the choice is usually limited by the size of the structure you want to study.  
CMA can offer a variety of probes with different membrane lengths from 1mm to 10 mm suitable for most experiments.

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<b>Question:</b>	Which membrane material is to be recommended? 怎样选择膜材料？
------------------	---

**Answer:** A membrane with low cut off purifies your sample by excluding large molecules. A membrane with high cut off recovers large substances such as peptides or smaller proteins. Note that some substances can also bind to the membrane material. To optimize your choice you should make an in vitro test with the substance you want to monitor. CMA can offer membranes made of PAES, polyarylethersulfone, 20,000 Daltons, cuprophane, 6,000 Daltons and PES polyetersulfone, 100,000 Daltons cut-off.

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**Answer:** There are of course many aspects which determines the choice. This is a brief overview, for more information see our catalog or product information on CMA web pages.

A stiff probe is suitable for a stereotaxic experiment on the brain while a flexible probe may be suited for microdialysis in a peripheral organ such as adipose tissue, muscle, liver or kidney.

CMA 12 the optimized probe for CNS use, ideal for chronic implantation

CMA 11 a thin stiff probe for discrete brain regions

CMA 7 an extremely small stiff probe for CNS studies in smaller animals such as mice

CMA 20 a soft non-metallic probe for peripheral tissues and blood vessels

CMA 30 a linear probe suitable for skin and other peripheral tissues

CMA 31 a linear probe suitable for skin and other peripheral tissues with 55k Da cut-off

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Question:	What perfusion flow rate should be used? 灌流液的流速是多少？
Answer:	<p>A high flow rate if you want to remove or introduce as many molecules as possible per time unit.</p> <p>A low flow rate should be used when you want to obtain a more concentrated dialysate (high recovery). Note that a low flow rate gives smaller volume. Consider also the volume needed for the analysis.</p> <p>如果你想在单位时间内去除或引入尽可能多的分子，应选用高流速。如果你想获得高浓度的透析液（高回收率），应选用低流速。但需要注意的是低流速收集的体积较小，因此在选用流速时还要考虑分析所需的体积。</p>
Question:	What time is needed to obtain steady state conditions? 什么时候才能达到稳态平衡？
Answer:	<p>The introduction of a probe into the tissue will always cause damage and the recovery of function will take a certain time period. An hour is often used to reach “baseline conditions”.</p> <p>当植入探针时会对组织造成一定的伤害，因此需要一段时间进行修复。通常灌流一个小时可达到“基线条件”。</p>
Question:	Are the probes reusable? 探针可以重复使用吗？
Answer:	<p>The preclinical probes can be used repeatedly if rinsed and stored in deionized water between experiments. However CMA can only guarantee single usage.</p> <p>实验中临床前研究使用的探针若清洗后储存在去离子水中，可以重复使用。但是，CMA 仅能保证可供一次使用。</p>

**Question:** When do I need the guides?  
什么时候需要探针导管？

**Answer:** Guides are used when you shall perform Microdialysis in the brain of a freely moving animal.

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**Question:** How do I handle FEP tubing and Tubing adapters?  
怎样连接微透析管路和管路接头？

When using a FEP tubing you shall cut the tubing with a sharp scalpel or similar to be sure the ends are open. After use rinse it with de ionized water to wash out the salts.

Tubing Adapters can be swelled in 70% ethanol for easy connections with tubing and syringes. They will shrink back in air again and ensure tight fit and zero internal volumes.

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Maybe you have not calibrated the pump to the syringe you are using.  
The syringe can be leaking or you are not using the syringe recommended by CMA.

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**Answer:** The probe may be leaking, the tubing may be blocked, something may be wrong with the pump, or the tubing adapters may be leaking.

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## Instruments

**Question:** What is the difference between the CMA 400 pump and the CMA 402 pump?  
CMA400 与 CMA402 区别有哪些？

**Answer:** CMA 400 has a pulse free flow from 1 nL/min - 1 mL/min, runs four syringes simultaneously. The pump is calibrated for various sizes of syringes. Apart from

microdialysis experiments, it can be used for microinjections of preset volumes that can be repeated in intervals.

CMA 402 has flow rates between 0,1µL - 20 µL/min, run by two syringes individually (1, 2.5 or 5 ml) . Start/stop and flow rate can be set independent for each syringe.

Both pumps can be used as both push and pull option.

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<b>Question:</b>	<p>What is the difference between CMA 142 and CMA 470 Fraction Collectors?</p> <p>CMA 142 和 CMA 470 微量收集器的区别有哪些？</p> <p>The CMA 142 collects fractions from 1 probe, 20 fractions or from 2 probes, 10 fractions. Fraction volumes from 1uL - 50 uL in open vials.</p> <p>CMA 470 is a refrigerated fraction collector, can collect 64 fractions from 1 up to 4 probes simultaneously (16 fractions/probe = 64 total). Fraction volumes from 1 uL - 2000 uL.</p> <p>Both opened and seal vials can be used.</p>
<b>Answer:</b>	<p>CMA 142 收集的样品量是 1 个探针，20 个收集瓶；或者 2 个探针，10 个收集瓶。透析液存储在开口小瓶，其体积范围是 1uL - 50 uL 。</p> <p>CMA 470 是低温微量收集器，可同时收集 1-4 个透析液，共收集 64 个样品。（每个探针收集 16 个样品，共 64 个）。透析液可存储在开口小瓶或密封小瓶中，其体积范围是 1 uL - 2000 uL</p>
<b>Question:</b>	<p>When shall the CMA 150 be used?</p> <p>CMA 150 什么时候使用？</p>
<b>Answer:</b>	<p>The CMA 150 keeps a stable temperature of the animal during an anesthetized microdialysis study. This ensures a better quality of your study.</p> <p>在微透析研究中 CMA 150 可维持麻醉动物的体温恒定，确保实验的准确性。</p>
<b>Question:</b>	<p>What is the use of a CMA 120 System for Freely Moving Animals?</p> <p>CMA 120 清醒活动系统的作用是什么？</p>
<b>Answer:</b>	<p>This system enables Microdialysis studies on conscious, small laboratory animals over a long period of time. It can be used in combination with any one of CMA's preclinical</p>

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