

# ORCHESTA™ AMBULATORY PUMP

Weighing only 180g, the OrchesTA model 500 pump (formerly Pegasus) is one of the most versatile pumps for ambulatory animal infusion. The peristaltic mechanism uses cost-efficient tube sets to deliver 0.1-100 ml/hr pulling from drug bags ranging from 50 to 300ml. The pumps are typically placed in the pocket of a Lomir jacket and the drug bag can be protected in a box.

Connect the model 500 pump with a cable to an OrchesTA pump transceiver and you can now control and monitor your infusion study remotely—limit human interaction with the animals, automate repetitive tasks and reduce the opportunities for programming and documentation errors.



Ideally suited for ambulatory infusion when placed in a jacket pocket. IV access is typically made through a subcutaneous port.



| Part No.    | Description   |
|-------------|---|
| OR-500-0001 | OrchesTA model 500 peristaltic pump (PEGA® LAB; 0.1-100ml/hr)       |
| OR-500-1000 | OrchesTA model 500 peristaltic pump tube (PEGA® tube 10255)         |
| OR-500-1205 | OrchesTA model 500 drug bag, PVC, 50mL (PEGA® bag 14050)            |
| OR-500-1210 | OrchesTA model 500 drug bag, PVC, 100mL (PEGA® bag 14100)           |
| OR-500-1215 | OrchesTA model 500 drug bag, PVC, 150mL (PEGA® bag 14150)           |
| OR-500-1230 | OrchesTA model 500 drug bag, PVC, 300mL (PEGA® bag 14300)           |
| OR-500-1405 | OrchesTA model 500 plastic box for 50ml drug bag (PEGA® box 10405)  |
| OR-500-1410 | OrchesTA model 500 plastic box for 100ml drug bag (PEGA® box 10410) |
| OR-500-1415 | OrchesTA model 500 plastic box for 150ml drug bag (PEGA® box 10415) |
| OR-500-1003 | OrchesTA PT500 wireless transceiver for model 500 pump              |
| OR-500-1080 | OrchesTA PCT200 PC transceiver, programmed for model 500 pumps      |

## SPECIFICATIONS

|              |  |
|--------------|--|
| Pump type    | Linear peristaltic                     |
| Flow rates   | 0.1-100ml/hr                           |
| Accuracy     | ±5%                                    |
| Pulse size   | 7µl                                    |
| Alarms       | Occlusion, bag empty, battery, others  |
| Alarm mute   | Yes                                    |
| Battery      | 2 x AA                                 |
| Battery life | ~1000ml with 3000mAh lithium batteries |
| Dimensions   | 8.7x6.4x3.3cm                          |
| Weight       | 180g with batteries                    |
| CE Mark      | Yes                                    |

NOTE: The OrchesTA model 500 version of this pump is not approved for human use.



Drug bags



Protective boxes attach to the bottom of the pump