

Minor Surgery Surgical Kit

SURGICAL INSTRUMENTS



ORDERING INFORMATION

Minor Surgery Surgical Kit

	Ref. No.
Minor Surgery Surgical Kit	72-8940
This kit consists of:	
Scalpel Handle No. 3	72-8350
Scalpel Blades No. 10, Sterile	72-8360
Eye Scissors, 11.5 cm Curved, Special Cut	72-8441
Blumenthal Bone Rongeurs, Curved, 15.5 cm	72-8906
Universal Clothing Scissors	72-8450
Operating Scissors, 14.5 cm Sharp/Blunt, Straight	72-8393
Kuehne Cover Glass Forceps, 10 cm, Angled	72-8676
Jeweller's Forceps 11 cm, No. 5 with Extra Delicate Points	72-8696
Probe 14.5 cm Double Ended, 1 mm Diameter	72-8912

Harvard Apparatus Surgical Instruments are forged and finished in a German ISO 9001 facility. These superior quality tools are highly resistant to corrosion and maintain their sharp cutting surfaces for a long time. You can purchase these surgical instruments individually or in a convenient kit.

The Minor Surgery Surgical Kit is ideally suited for basic procedures and minor surgery.

Major Deluxe Surgical Kit

SURGICAL INSTRUMENTS



Harvard Apparatus Surgical Instruments are forged and finished in a German ISO 9001 facility. These superior quality tools are highly resistant to corrosion and maintain their sharp cutting surfaces for a long time. You can purchase these surgical instruments individually or in a convenient kit.

The Major Deluxe Surgical Kit is ideally suited for more extensive procedures and major laboratory animal surgeries.

ORDERING INFORMATION

Major Deluxe Surgical Kit	
	Ref. No.
Major Deluxe Surgical Kit	72-8941
This kit consists of:	
Scalpel Handle No. 3	72-8350
Scalpel Blades No. 10, Sterile	72-8360
Scalpel Blades No. 11, Sterile	72-8362
Vannas Eye Scissors, Spring Action Model Tübingen, 8.5 cm, Straight	72-8508
Olsen-Hegar Needle Holders with Scissors with Tungsten Carbide Inserts, 14.0 cm	72-8899
Eye Scissors, 11.5 cm, Straight, Special Cut	72-8440
Jeweller's Forceps, 11.0 cm, No. 5, with Extra Delicate Points	72-8696
Graefe Iris Forceps, Serrated, 7.0 cm, Curved, Points 0.5 mm	72-8598
UltraEdge Operating Scissors, 14.5 cm Sharp/Blunt, Straight	72-8394
Universal Clothing Scissors	72-8450
Friedman Bone Rongeurs, Delicate	72-8905
Adson Tissue Forceps, 12.0 cm, 1x2 Teeth	72-8547
Barraquer (Colibri) Eye Specula, Large	72-8919
Probes, 14.5 cm, Double Ended, 1 mm, Diameter	72-8912

Rat Surgical Kit

SURGICAL INSTRUMENTS



Harvard Apparatus Surgical Instruments are forged and finished in a German ISO 9001 facility. These superior quality tools are highly resistant to corrosion and maintain their sharp cutting surfaces for long time. You can purchase these surgical instruments individually or in a convenient kit.

The Rat Surgical Kit is an ideal starting point for surgical procedures commonly performed on laboratory rats.

ORDERING INFORMATION

Rat Surgical Kit	Ref. No.
Rat Surgical Kit	72-8942
This kit consists of:	
Disposable Scalpels with plastic handle No. 11, sterile box of 10	72-8372
Jeweler's Forceps 11.0 cm No. 5, with Extra Delicate Points	72-8696
Vannas Eye Scissors 8.5 cm Straight, Spring Action	72-8503
Graefe Iris Forceps Serrated 10.0 cm Curved Points, 0.7 mm	72-8604
Eye Scissors 11.5 cm Straight, Special Cut	72-8440
Barraquer (Colibri) Eye Specula, Large	72-8919
Olsen-Hegar Needle Holder with Scissors, 11.5 cm	72-8844
Probes 14.5 cm Double Ended, 1 mm Diameter	72-8912

Mouse Surgical Kit

SURGICAL INSTRUMENTS



Harvard Apparatus Surgical Instruments are forged and finished in a German ISO 9001 facility. These superior quality tools are highly resistant to corrosion and maintain their sharp cutting surfaces for a long time. You can purchase these surgical instruments individually or in a convenient kit.

The Mouse Surgical Kit is ideal starting point for surgical procedures performed on laboratory mice.

ORDERING INFORMATION

Mouse Surgical Kit	Ref. No.
Mouse Surgical Kit	72-8943
This kit consists of:	
Disposable Scalpels with Plastic Handle No. 11, Sterile, Box of 10	72-8372
Jeweller's Forceps No. 5, 11.0 cm with Extra Delicate Points	72-8696
Vannas Eye Scissors, 8.5 cm Straight Spring Action	72-8503
Graefe Iris Forceps Serrated, 10.0 cm, Curved Points 0.7 mm	72-8604
Eye Scissors 10.5 cm, Straight	72-8426
Vannas Eye Scissors, 8.5 cm Barraquer (Colibri) Eye Specula, Small	72-8918
Derf Needle Holders, 12.0 cm	72-8826
Probe, 14.5 cm Double Ended, 1 mm Diameter	72-8912