


ZensorSPE


網版印刷碳電極

Screen-Printed Electrodes (Carbon)




Performance

GCE



TE100



Condition

Solution : 3mM Ferricyanide, 0.1M KCl

Specification

Type	Appearance	Model	Description (Diameter/Area)
Single Electrode		SE100	WE : 5mm / 0.196cm ²
		SE101	WE : 3mm / 0.071cm ²
		SE102	WE : 1.5mm / 0.018cm ²
Three Electrode		TE100	WE : 3mm / 0.071cm ² CE : Carbon (0.050cm ²) RE : Ag (0.010cm ²)
		TE200	WE : 2x1mm / 0.020cm ² CE : Carbon (0.040cm ²) RE : Ag (0.020cm ²)

